**READ ALL INSTRUCTIONS
BEFORE HANDLING DEVICE**

Device Description

ENGLISH

HemaXis

HemaXis DB 5µL and 10µL Capillary Blood Collection Procedure

- HemaXis DB is intended to collect, store and transport capillary blood samples for analysis in a certified laboratory. HemaXis DB is recommended for use under professional supervision in a medical environment. Single use only
- To avoid contamination risk, use device directly after removing from package.
- The HemaXis DB device may be stored at room temperature.


EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- Capillary Blood collection supplies: sterile retractable lancets, alcohol swabs, sterile gauze and adhesive bandages/plasters.
- Personal protective equipment if assisting patient to collect blood: Gloves and other personal protective equipment as required to provide protection from exposure to blood borne pathogens.


OPTIONAL MATERIALS NOT PROVIDED

- Mailing envelope: for transport of specimens.
- Impermeable bag and desiccant: for transport of specimens.

Watch instructional video at www.hemaxis.com


10µL channels are curved and 5µL channels are straight.

- Blood samples can only be applied and collected through the Channel Entrances along the edge of the collection chip.
- Open kit and place on a flat surface before using.

1 Open kit and draw blood from the fingertip

- a. Sterilize index or ring finger with alcohol swap.
- b. Prick finger with a sterile lancet.
- c. Wipe off first blood drop with tissue (not alcohol swap) then gently press finger to produce another drop of blood.

2 Fill ONE collection channel with blood

Small blood meniscus forms at channel outlet hole

Only fill ONE CHANNEL with blood then proceed to next step to transfer. DO NOT fill more than one channel at a time.

- a. Position fingertip just below the channel entrance and make contact with the blood drop to the channel entrance until the channel is completely filled and meniscus appears at the channel outlet.
- b. Wipe excess blood from finger using sterile gauze or tissue

3 Transfer blood to card - Repeat 3 times

- a. Firmly close the device for **10 seconds** to transfer blood sample to the filter card. The closure flap doesn't need to be in locked position.
- b. A blood spot will be seen forming on the filter card through one of the round control windows.
- c. Repeat steps 2 and 3 until four blood spots are generated.

- Dry samples for 60 min. then close device and lock closure.
- Store or ship device according to your procedures.

Troubleshooting Tips

✗ Not enough blood volume forms at the fingertip to fill the collection channel

- ✓ Wash patient's hands in very warm water.
- ✓ Shake and massage fingers of selected hand to promote blood flow.
- ✓ Use a lancet with a longer needle. We recommend using a 1.8 mm needle.
- ✓ Prick finger again with a new lancet

✗ Collection channel does not fill completely

- ✓ Reapply finger to the inlet 2-3 times to completely fill and ensure meniscus forms at the channel outlet.
- ✓ Hold the device with the Collection Panel angled upwards to allow gravity to help pull blood through the channel
- ✓ Prick finger again with a new lancet to generate more blood

✗ Blood does not transfer to filter paper card

- ✓ Open device, and reapply finger to the inlet 2-3 times to completely fill and ensure meniscus forms
- ✓ Blood may have coagulated in the channel. The transfer must be done within 30 seconds following the filling of a channel.
- ✓ The device must remain closed for a minimum of 10 seconds to enable the complete transfer of the blood to the filter card.

✗ Erratic blood spotting observed on the filter paper card

- ✓ Too much blood was applied to the channel inlet. Do not press fingertip against the inlet so that skin touches the edge.
- ✓ The blood spots in the correct location can still be used as long as there is more than 2 mm between the main spot and another spot.
- ✓ Do not attempt to refill a used collection channel.

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