**READ ALL INSTRUCTIONS
BEFORE HANDLING DEVICE**

Device Description

ENGLISH

HemaXis

HemaXis DB 5µL and 10µL Capillary Blood Collection Procedure

- HemaXis DB is intended to collect, store and transport capillary blood samples for analysis in a certified laboratory. HemaXis DB is recommended for use under professional supervision in a medical environment. Single use only.
- To avoid contamination risk, use device directly after removing from package.
- The HemaXis DB device may be stored at room temperature.


EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- Capillary Blood collection supplies: sterile retractable lancets, alcohol swabs, sterile gauze and adhesive bandages/plasters.
- Personal protective equipment if assisting patient to collect blood: Gloves and other personal protective equipment as required to provide protection from exposure to blood borne pathogens.


OPTIONAL MATERIALS NOT PROVIDED

- Mailing envelope: for transport of specimens.
- Impermeable bag and desiccant: for transport of specimens.

Watch instructional video at www.hemaxis.com/products/documentation


10µL channels are curved and 5µL channels are straight.

- Blood samples can only be applied and collected through the Channel Entrances along the edge of the collection chip.
- Open kit and place on a flat surface before using.

1 Generate blood drop at the fingertip

- Sterilize index or ring finger with alcohol swab
- Prick fingertip with a sterile lancet
- Wipe off first blood drop with tissue (not alcohol swab)
- Gently press finger to produce another drop of blood

2 Fill ONE collection channel with blood


Small blood meniscus forms at channel outlet hole

Only fill ONE CHANNEL at a time

- Open the device and bend the closure flap completely backward
- Position fingertip just below the channel entrance and make contact with the blood drop to the channel entrance until the channel is completely filled. **Always position the device so that the channel outlet is lower than the channel inlet**
- A meniscus must appear at the channel outlet. You may need to make contact 2-3 times to ensure complete filling of the channel
- Wipe excess blood from finger using sterile gauze or tissue

3 Transfer blood to card - Repeat 3 times

Control window

- Close the device for **10 seconds** to transfer blood sample to the filter card. The closure flap doesn't need to be in locked position
- A blood spot will be seen forming on the filter card through one of the round control windows on the filter card panel
- Repeat steps 2 and 3 until four blood spots are generated

- Let the samples dry for 60 min. then close the device and lock closure flap
- Store or ship device according to your procedures

Troubleshooting Tips

✗ Not enough blood volume forms at the fingertip to fill the collection channel

- ✓ Wash patient's hands in very warm water
- ✓ Shake and massage fingers to activate blood flow
- ✓ Use a lancet with a longer needle. We recommend using a 1.8 mm needle
- ✓ Prick finger again with a new lancet

✗ Collection channel does not fill completely

- ✓ Make sure the channel outlet is maintain lower than the channel inlet, to allow gravity to help pull blood through the channel
- ✓ Reapply finger to the inlet 2-3 times to completely fill and ensure meniscus forms at the channel outlet
- ✓ Prick finger again with a new lancet to generate more blood

✗ Blood does not transfer to filter paper card

- ✓ Open device, and reapply finger to the inlet 2-3 times to completely fill and ensure meniscus forms
- ✓ The device must remained closed for a minimum of 10 seconds to enable the complete transfer of the blood to the filter card
- ✓ Blood may have coagulated in the channel. The transfer must be done within 30 seconds following the filling of a channel

✗ Erratic blood spotting observed on the filter paper card

- ✓ Too much blood was applied to the channel inlet. Do not press fingertip against the inlet so that skin touches the edge
- ✓ The blood spots in the correct location can still be used as long as there is more than 2 mm between the main spot and another spot
- ✓ Do not attempt to refill a used collection channel

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